

## Product Datasheet

# Canine Thyroid Stimulating Hormone (TSH) ELISA Kit (orb1203945)

<b>Description</b>	Canine Thyroid Stimulating Hormone (TSH) ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	46.8-3000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9XSW8</a>
<b>Sensitivity</b>	11.25pg/mL
<b>Expiration Date</b>	Please enquire.