

## Product Datasheet

# Mouse Transient Receptor Potential Cation Channel Subfamily V, Member 2 (TRPV2) ELISA Kit (orb1203946)

<b>Description</b>	Mouse Transient Receptor Potential Cation Channel Subfamily V, Member 2 (TRPV2) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.156-10ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9WTR1</a>
<b>Sensitivity</b>	0.049ng/mL
<b>Entrez</b>	<a href="#">22368</a>
<b>Expiration Date</b>	Please enquire.