

## Product Datasheet

# Canine Surfactant Associated Protein A (SPA) ELISA Kit (orb1203999)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Canine Surfactant Associated Protein A (SPA) ELISA Kit                           |
| <b>Reactivity</b>        | Canine   |
| <b>Range</b>             | 1.87-120ng/mL  |
| <b>Note</b>              | For research use only  |
| <b>Application notes</b> | sandwich   |
| <b>Sample Types</b>      | serum, plasma, lung lavage fluid, tissue homogenates or other biological fluids. |
| <b>Uniprot ID</b>        | <a href="#">Q6RXL1</a>   |
| <b>Sensitivity</b>       | 0.62ng/mL  |
| <b>Entrez</b>            | <a href="#">407213</a>   |
| <b>Expiration Date</b>   | Please enquire.  |