

Product Datasheet

Mouse Secreted Frizzled Related Protein 4 (SFRP4) ELISA Kit (orb1204026)

Description	Mouse Secreted Frizzled Related Protein 4 (SFRP4) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q9Z1N6
Sensitivity	0.065ng/mL
Entrez	20379
Expiration Date	Please enquire.