

Product Datasheet

Mouse Regenerating Islet Derived Protein 3 Beta (REG3b) ELISA Kit (orb1204057)

Description	Mouse Regenerating Islet Derived Protein 3 Beta (REG3b) ELISA Kit
Reactivity	Mouse
Range	9.37-600pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	P35230
Sensitivity	4.1pg/mL
Entrez	18489
Expiration Date	Please enquire.