

Product Datasheet

Mouse Regenerating Islet Derived Protein 3 Alpha (REG3a) ELISA Kit (orb1204058)

Description	Mouse Regenerating Islet Derived Protein 3 Alpha (REG3a) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	O09037
Sensitivity	6.3pg/mL
Entrez	19694
Expiration Date	Please enquire.