

## Product Datasheet

# Guinea pig Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit (orb1204060)

<b>Description</b>	Guinea pig Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit
<b>Reactivity</b>	Guinea pig
<b>Range</b>	0.156-10ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">P97272</a>
<b>Sensitivity</b>	0.059ng/mL
<b>Entrez</b>	<a href="#">100135495</a>
<b>Expiration Date</b>	Please enquire.