

Product Datasheet

Canine Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit (orb1204061)

Description	Canine Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit
Reactivity	Canine
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8HYS0
Sensitivity	7.8pg/mL
Entrez	403522
Expiration Date	Please enquire.