

## Product Datasheet

# Mouse Peroxiredoxin 4 (PRDX4) ELISA Kit (orb1204073)

<b>Description</b>	Mouse Peroxiredoxin 4 (PRDX4) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.312-20ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<b>008807</b>
<b>Sensitivity</b>	0.108ng/mL
<b>Entrez</b>	<b>53381</b>
<b>Expiration Date</b>	Please enquire.