

Product Datasheet

Mouse Peripheral Myelin Protein 22 (PMP22) ELISA Kit (orb1204084)

Description	Mouse Peripheral Myelin Protein 22 (PMP22) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P16646
Sensitivity	0.058ng/mL
Entrez	18858
Expiration Date	Please enquire.