

## Product Datasheet

# Mouse Protein Disulfide Isomerase A4 (PDIA4) ELISA Kit (orb1204096)

<b>Description</b>	Mouse Protein Disulfide Isomerase A4 (PDIA4) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.312-20ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<b>P08003</b>
<b>Sensitivity</b>	0.112ng/mL
<b>Entrez</b>	<b>12304</b>
<b>Expiration Date</b>	Please enquire.