

Product Datasheet

Mouse Prenylcysteine Oxidase 1 (PCYOX1) ELISA Kit (orb1204100)

Description	Mouse Prenylcysteine Oxidase 1 (PCYOX1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9CQF9
Sensitivity	0.072ng/mL
Entrez	66881
Expiration Date	Please enquire.