

Product Datasheet

Mouse Neuropeptide Y Receptor Y2 (NPY2R) ELISA Kit (orb1204123)

Description	Mouse Neuropeptide Y Receptor Y2 (NPY2R) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P97295
Sensitivity	0.053ng/mL
Entrez	18167
Expiration Date	Please enquire.