

Product Datasheet

Canine Matrix Metalloproteinase 3 (MMP3) ELISA Kit (orb1204158)

Description	Canine Matrix Metalloproteinase 3 (MMP3) ELISA Kit
Reactivity	Canine
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q6Y4Q5
Sensitivity	0.063ng/mL
Entrez	403445
Expiration Date	Please enquire.