

Product Datasheet

Canine Macrophage Inflammatory Protein 3 Beta (MIP3b) ELISA Kit (orb1204160)

Description	Canine Macrophage Inflammatory Protein 3 Beta (MIP3b) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q68AY7
Sensitivity	6.3pg/mL
Entrez	448793
Expiration Date	Please enquire.