

Product Datasheet

Canine Macrophage Inflammatory Protein 1 Alpha (MIP1a) ELISA Kit (orb1204162)

Description	Canine Macrophage Inflammatory Protein 1 Alpha (MIP1a) ELISA Kit
Reactivity	Canine
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q68A92
Sensitivity	28pg/mL
Entrez	448787
Expiration Date	Please enquire.