

## Product Datasheet

# Canine Monocyte Chemotactic Protein 2 (MCP2) ELISA Kit (orb1204169)

<b>Description</b>	Canine Monocyte Chemotactic Protein 2 (MCP2) ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	31.2-2000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q68AY9</a>
<b>Sensitivity</b>	9.7pg/mL
<b>Entrez</b>	<a href="#">448792</a>
<b>Expiration Date</b>	Please enquire.