

Product Datasheet

Canine Leukemia Inhibitory Factor (LIF) ELISA Kit (orb1204190)

Description	Canine Leukemia Inhibitory Factor (LIF) ELISA Kit
Reactivity	Canine
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	J9NRL6
Sensitivity	14.7pg/mL
Entrez	403449
Expiration Date	Please enquire.