

Product Datasheet

Mouse Insulin Promoter Factor 1 (IPF) ELISA Kit (orb1204219)

Description	Mouse Insulin Promoter Factor 1 (IPF) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P52946
Sensitivity	0.053ng/mL
Entrez	18609
Expiration Date	Please enquire.