

Product Datasheet

Mouse Insulin Like Protein 6 (INSL6) ELISA Kit (orb1204220)

Description	Mouse Insulin Like Protein 6 (INSL6) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9QY05
Sensitivity	30pg/mL
Entrez	27356
Expiration Date	Please enquire.