

Product Datasheet

Mouse Heparan Sulfate Proteoglycan 2 (HSPG2) ELISA Kit (orb1204234)

Description	Mouse Heparan Sulfate Proteoglycan 2 (HSPG2) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q05793
Sensitivity	32.4pg/mL
Entrez	15530
Expiration Date	Please enquire.