

Product Datasheet

Mouse Hyaluronan Synthase 2 (HAS2) ELISA Kit (orb1204253)

Description	Mouse Hyaluronan Synthase 2 (HAS2) ELISA Kit
Reactivity	Mouse
Range	78-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P70312
Sensitivity	26pg/mL
Entrez	15117
Expiration Date	Please enquire.