

Product Datasheet

Mouse Glutamate Receptor, Ionotropic, AMPA 2 (GRIA2) ELISA Kit (orb1204264)

Description	Mouse Glutamate Receptor, Ionotropic, AMPA 2 (GRIA2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P23819
Sensitivity	0.058ng/mL
Entrez	14800
Expiration Date	Please enquire.