

## Product Datasheet

### Mouse Fibulin 3 (FBLN3) ELISA Kit (orb1204306)

<b>Description</b>	Mouse Fibulin 3 (FBLN3) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	15.6-1000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, urine or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q8BPB5</a>
<b>Sensitivity</b>	6.6pg/mL
<b>Entrez</b>	<a href="#">216616</a>
<b>Expiration Date</b>	Please enquire.