

Product Datasheet

Guinea pig Eosinophil Chemotactic Factor (ECF) ELISA Kit (orb1204330)

Description	Guinea pig Eosinophil Chemotactic Factor (ECF) ELISA Kit
Reactivity	Guinea pig
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	G8XRI2
Sensitivity	6.2pg/mL
Entrez	100038010
Expiration Date	Please enquire.