

Product Datasheet

Mouse Epstein Barr Virus Induced Protein 3 (EBI3) ELISA Kit (orb1204331)

Description	Mouse Epstein Barr Virus Induced Protein 3 (EBI3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	O35228
Sensitivity	0.054ng/mL
Entrez	50498
Expiration Date	Please enquire.