

## Product Datasheet

# Mouse Cytochrome P450 27B1 (CYP27B1) ELISA Kit (orb1204364)

<b>Description</b>	Mouse Cytochrome P450 27B1 (CYP27B1) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.625-40ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates, cell lysates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">O35084</a>
<b>Sensitivity</b>	0.214ng/mL
<b>Entrez</b>	<a href="#">13115</a>
<b>Expiration Date</b>	Please enquire.