

Product Datasheet

Mouse Calcium Channel, Voltage Dependent, L-Type, Alpha 1S Subunit (CACNa1S) ELISA Kit (orb1204436)

Description	Mouse Calcium Channel, Voltage Dependent, L-Type, Alpha 1S Subunit (CACNa1S) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q02789
Sensitivity	0.101ng/mL
Entrez	12292
Expiration Date	Please enquire.