

## Product Datasheet

# Mouse ATPase, H<sup>+</sup> Transporting, Lysosomal Accessory Protein 2 (ATP6AP2) ELISA Kit (orb1204454)

<b>Description</b>	Mouse ATPase, H <sup>+</sup> Transporting, Lysosomal Accessory Protein 2 (ATP6AP2) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.156-10ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9CYN9</a>
<b>Sensitivity</b>	0.055ng/mL
<b>Entrez</b>	<a href="#">70495</a>
<b>Expiration Date</b>	Please enquire.