

Product Datasheet

Monkey Interleukin 4 (IL4) ELISA Kit (orb1204516)

Description	Monkey Interleukin 4 (IL4) ELISA Kit
Reactivity	Monkey
Range	3.12-200pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P51492
Sensitivity	1.26pg/mL
Entrez	574281
Expiration Date	Please enquire.