

Product Datasheet

Human Cysteine-Rich Secretory Protein LCCL Domain-Containing 2 (CRISPLD2) ELISA Kit (orb1204534)

Description	Human Cysteine-Rich Secretory Protein LCCL Domain-Containing 2 (CRISPLD2) ELISA Kit
Reactivity	Human
Range	3.12-200ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9H0B8
Sensitivity	1.12ng/mL
Entrez	83716
Expiration Date	Please enquire.