

Product Datasheet

Human 5-Hydroxytryptamine Receptor 7 (HTR7) ELISA Kit (orb1204552)

Description	Human 5-Hydroxytryptamine Receptor 7 (HTR7) ELISA Kit
Reactivity	Human
Range	78-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P34969
Sensitivity	34.2pg/mL
Entrez	3363
Expiration Date	Please enquire.