

## Product Datasheet

# Human WNT1 Inducible Signaling Pathway Protein 2 (WISP2) ELISA Kit (orb1204583)

<b>Description</b>	Human WNT1 Inducible Signaling Pathway Protein 2 (WISP2) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312-20ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">O76076</a>
<b>Sensitivity</b>	0.133ng/mL
<b>Entrez</b>	<a href="#">8839</a>
<b>Expiration Date</b>	Please enquire.