

Product Datasheet

Human VEGF Co Regulated Chemokine 1 (VCC1) ELISA Kit (orb1204598)

Description	Human VEGF Co Regulated Chemokine 1 (VCC1) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q6UXB2
Sensitivity	4.7pg/mL
Entrez	284340
Expiration Date	Please enquire.