

Product Datasheet

Human Toll Like Receptor Adaptor Molecule 1 (TICAM1) ELISA Kit (orb1204692)

Description	Human Toll Like Receptor Adaptor Molecule 1 (TICAM1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q8IUC6
Sensitivity	0.124ng/mL
Entrez	148022
Expiration Date	Please enquire.