

Product Datasheet

Human Six Transmembrane Epithelial Antigen Of The Prostate 2 (STEAP2) ELISA Kit (orb1204745)

Description	Human Six Transmembrane Epithelial Antigen Of The Prostate 2 (STEAP2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q8NFT2
Sensitivity	0.062ng/mL
Entrez	261729
Expiration Date	Please enquire.