

Product Datasheet

Human Sclerostin Domain Containing Protein 1 (SOSTDC1) ELISA Kit (orb1204772)

Description	Human Sclerostin Domain Containing Protein 1 (SOSTDC1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q6X4U4
Sensitivity	0.059ng/mL
Entrez	25928
Expiration Date	Please enquire.