

Product Datasheet

Human Mothers Against Decapentaplegic Homolog 3 (Smad3) ELISA Kit (orb1204787)

Description	Human Mothers Against Decapentaplegic Homolog 3 (Smad3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P84022
Sensitivity	0.066ng/mL
Entrez	4088
Expiration Date	Please enquire.