

Product Datasheet

Human Mothers Against Decapentaplegic Homolog 2 (Smad2) ELISA Kit (orb1204788)

Description	Human Mothers Against Decapentaplegic Homolog 2 (Smad2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q15796
Sensitivity	0.066ng/mL
Entrez	4087
Expiration Date	Please enquire.