

Product Datasheet

Human Solute Carrier Family 39, Member 6 (SLC39A6) ELISA Kit (orb1204793)

Description	Human Solute Carrier Family 39, Member 6 (SLC39A6) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q13433
Sensitivity	0.114ng/mL
Entrez	25800
Expiration Date	Please enquire.