

Product Datasheet

Human Sirtuin 6 (SIRT6) ELISA Kit (orb1204808)

Description	Human Sirtuin 6 (SIRT6) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8N6T7
Sensitivity	0.124ng/mL
Entrez	51548
Expiration Date	Please enquire.