

Product Datasheet

Human Amiloride Sensitive Sodium Channel Subunit Alpha (SCNN1a) ELISA Kit (orb1204848)

Description	Human Amiloride Sensitive Sodium Channel Subunit Alpha (SCNN1a) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P37088
Sensitivity	0.067ng/mL
Entrez	6337
Expiration Date	Please enquire.