

Product Datasheet

Human SATB Homeobox Protein 1 (SATB1) ELISA Kit (orb1204858)

| | |
|--------------------------|---|
| Description | Human SATB Homeobox Protein 1 (SATB1) ELISA Kit |
| Reactivity | Human |
| Range | 0.312-20ng/mL |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | tissue homogenates or other biological fluids. |
| Uniprot ID | Q01826 |
| Sensitivity | 0.102ng/mL |
| Entrez | 6304 |
| Expiration Date | Please enquire. |