

Product Datasheet

Human Relaxin/Insulin Like Family Peptide Receptor 3 (RXFP3) ELISA Kit (orb1204871)

Description	Human Relaxin/Insulin Like Family Peptide Receptor 3 (RXFP3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9NSD7
Sensitivity	0.057ng/mL
Entrez	51289
Expiration Date	Please enquire.