

Product Datasheet

Human Proteasome Subunit Alpha Type 1 (PSMa1) ELISA Kit (orb1204944)

Description	Human Proteasome Subunit Alpha Type 1 (PSMa1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P25786
Sensitivity	0.067ng/mL
Entrez	5682
Expiration Date	Please enquire.