

Product Datasheet

Human Phosphatidylinositol Binding Clathrin Assembly Protein (PICALM) ELISA Kit (orb1205019)

Description	Human Phosphatidylinositol Binding Clathrin Assembly Protein (PICALM) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q13492
Sensitivity	0.126ng/mL
Entrez	8301
Expiration Date	Please enquire.