

Product Datasheet

Human Protein Disulfide Isomerase A4 (PDIA4) ELISA Kit (orb1205035)

Description	Human Protein Disulfide Isomerase A4 (PDIA4) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P13667
Sensitivity	0.66ng/mL
Entrez	9601
Expiration Date	Please enquire.