

Product Datasheet

Human Purinergic Receptor P2Y, G Protein Coupled 12 (P2RY12) ELISA Kit (orb1205076)

Description	Human Purinergic Receptor P2Y, G Protein Coupled 12 (P2RY12) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9H244
Sensitivity	0.063ng/mL
Entrez	64805
Expiration Date	Please enquire.