

Product Datasheet

Human Neuropeptide FF (NPFF) ELISA Kit (orb1205131)

Description	Human Neuropeptide FF (NPFF) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	015130
Sensitivity	4.54pg/mL
Entrez	8620
Expiration Date	Please enquire.