

Product Datasheet

Human Myosin Heavy Chain 3, Skeletal Muscle, Embryonic (MYH3) ELISA Kit (orb1205190)

Description	Human Myosin Heavy Chain 3, Skeletal Muscle, Embryonic (MYH3) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P11055
Sensitivity	0.126ng/mL
Entrez	4621
Expiration Date	Please enquire.